

Hemostatic Pellets

Keep your cases moving!

Aid in gingival retraction and hemostasis.



- Soft pellets are activated by normal tissue moisture
- For rapid activation, dip in water just prior to application
- Packaged in a spring-loaded, self-dispensing weighted plastic container

Epidri Pellets

- Cotton Pellets with 1.9mg Racemic Epinephrine HCl per pellet
- Item Number 12-100
- Color coded yellow




Racellet Pellets

- Size #2 1.5mg Racemic Epinephrine HCl per pellet
- Item Number 12-150
- Size #3 0.55mg Racemic Epinephrine HCl per pellet
- Item Number 12-200
- Color coded green



Retrax® Pellets

- Cotton Pellets with 15% Aluminum Sulfate per pellet
- Item Number 12-400
- Color coded orange




Packing Instruments



MOST POPULAR

R-55 - Serrated



R-55A - Serrated



R-50 - Plain



Guyer 7A - Serrated




R-101 - Plain




R-11 - Serrated



Retraction Cord Comparison Chart



Knitted



Braided



Twisted

37 varieties of retraction cord.



Global Leader
in Tissue Management Products



Please, contact
Pascal International or
your dealer for more
information.

ph **800.426.8051**
www.pascalinternational.net

Knitted Retraction Cord

Compressible interlocking chains transport greater amounts of hemostatic agent to the sulcus, offering more predictable, effective hemostasis and tissue management.

Available in:

KNITTRAX™ - Unmedicated

Braided Retraction Cord

Braided cotton cord strands will not separate, and packing instruments will not penetrate and seize the cord strands.

Available in:

SIL-TRAX PLUS® - Epinephrine/Zinc Phenosulfonate

SIL-TRAX EPI® - Epinephrine

SIL-TRAX AS® - Aluminum Sulfate

SIL-TRAX PLAIN® - Unmedicated

Twisted Retraction Cord

Soft, absorbent cord of long fiber cotton. Tightly wound to prevent unraveling and fraying during packing.

Available in:

RACORD TWO® - Epinephrine/Zinc Phenosulfonate

RACORD® - Epinephrine

PASCORD® - Aluminum Sulfate

RETRAX® - Unmedicated

Cord comparison

Pascal Cord Type Knitted Non-Medicated	GingiPak Equivalent	Ultradent Equivalent	Kerr/Dux Equivalent
Knittrax 07-569 Size 000 Non- Impregnated	none	Ultrapak Size 000 Non- Impregnated	GingiKNIT Size 000n Non- Impregnated
Knittrax 07-570 Size 00 Non- Impregnated	none	Ultrapak Size 00 Non- Impregnated	GingiKNIT Size 00n Non- Impregnated
Knittrax 07-575 Size 0 Non- Impregnated	none	Ultrapak Size 0 Non- Impregnated	GingiKNIT Size 0n Non- Impregnated
Knittrax 07-580 Size 1 Non- Impregnated	none	Ultrapak Size 1 Non- Impregnated	GingiKNIT Size 1n Non- Impregnated
Knittrax 07-585 Size 2 Non- Impregnated	none	Ultrapak Size 2 Non- Impregnated	GingiKNIT Size 2n Non- Impregnated

Pascal sample cords are unmedicated. Actual cords are color coded by size and/or medication.

Pascal Cord Type Braided w/ Epinephrine	GingiPak Equivalent	Kerr/Dux Equivalent
Siltrax EPI 07-197 Size 7/00 Epinephrine	Gingi-Pak Z-Twist Max 00	none
Siltrax EPI 07-200 Size 8/0 Epinephrine	Gingi-Pak Z-Twist Max 1	none
Siltrax EPI 07-210 Size 9/1 Epinephrine	Gingi-Pak Z-Twist Max 2	none
Siltrax EPI 07-220 Size 10/2 Epinephrine	Gingi-Pak Z-Twist Max 3	none

Pascal Cord Type Twisted w/ Epinephrine	GingiPak Equivalent	Kerr/Dux Equivalent
Racord 07-697 Size 7/00 Epinephrine	none	none
Racord 07-700 Size 8/0 Epinephrine	Gingi-Pak Max Size 1 Epinephrine	none
Racord 07-710 Size 9/1 Epinephrine	Gingi-Pak Max Size 2 Epinephrine	none
Racord 07-720 Size 10/2 Epinephrine	Gingi-Pak Max Size 3 Epinephrine	none

Pascal Cord Type Braided Aluminum Sulfate	GingiPak Equivalent	Kerr/Dux Equivalent
Siltax AS 07-097 Size 7/00 Alum. Sulfate	Gingi-Aid Z-Twist 00	GingiBraid + or Shortcut #0a Alum. Pot. Sulfate
Siltax AS 07-100 Size 8/0 Alum. Sulfate	Gingi-Aid Z-Twist 1	GingiBraid + or Shortcut #1a Alum. Pot. Sulfate
Siltax AS 07-110 Size 9/1 Alum. Sulfate	Gingi-Aid Z-Twist 2	GingiBraid + or Shortcut #2a Alum. Pot. Sulfate
Siltax AS 07-120 Size 10/2 Alum. Sulfate	Gingi-Aid Z-Twist 3	GingiBraid + or Shortcut #3a Alum. Pot. Sulfate

Pascal Cord Type Twisted Epi/Astringent	GingiPak Equivalent	Kerr/Dux Equivalent
Racord II 07-797 Size 7/00 Epi/Astringent	none	none
Racord II 07-800 Size 8/0 Epi/Astringent	none	GingiCORD #1 Epi/Alum
Racord II 07-810 Size 9/1 Epi/Astringent	none	GingiCORD #2 Epi/Alum
Racord II 07-820 Size 10/2 Epi/Astringent	none	GingiCORD #3 Epi/Alum

Pascal Cord Type Braided Epi/Astringent	GingiPak Equivalent	Kerr/Dux Equivalent
Siltrax Plus 07-297 Size 7/00 Epi/Astringent	none	GingiBraid + or Shortcut #0e Epi/Alum
Siltrax Plus 07-300 Size 8/0 Epi/Astringent	none	GingiBraid + or Shortcut #1e Epi/Alum
Siltrax Plus 07-310 Size 9/1 Epi/Astringent	none	GingiBraid + or Shortcut #2e Epi/Alum
Siltrax Plus 07-320 Size 10/2 Epi/Astringent	none	GingiBraid + or Shortcut #3e Epi/Alum

Pascal Cord Type Braided Non-Medicated	GingiPak Equivalent	Kerr/Dux Equivalent
Siltrax Plain 07-397 Size 7/00 Non-Impregnated	Gingi-Plain Max Z-Twist 00	GingiBraid + Shortcut #0n Non-Impregnated
Siltrax Plain 07-400 Size 8/0 Non-Impregnated	Gingi-Plain Max Z-Twist 1	GingiBraid + Shortcut #1n Non-Impregnated
Siltrax Plain 07-410 Size 9/1 Non-Impregnated	Gingi-Plain Max Z-Twist 2	GingiBraid + Shortcut #2n Non-Impregnated
Siltrax Plain 07-420 Size 10/2 Non-Impregnated	Gingi-Plain Max Z-Twist 3	GingiBraid + Shortcut #3n Non-Impregnated

Pascal Cord Type Twisted w/Alum. Sulfate	GingiPak Equivalent	Kerr/Dux Equivalent
Pascord 07-597 Size 7/00 Alum. Sulfate	none	none
Pascord 07-600 Size 8/0 Alum. Sulfate	none	none
Pascord 07-610 Size 9/1 Alum. Sulfate	none	none
Pascord 07-620 Size 10/2 Alum. Sulfate	none	none

Pascal Cord Type Twisted Non-Medicated	GingiPak Equivalent	Kerr/Dux Equivalent
Retrax 07-897 Size 7/00 Non-Impregnated	none	none
Retrax 07-900 Size 8/0 Non-Impregnated	Gingi-Plain #1 Max Soft Twist	none
Retrax 07-910 Size 9/1 Non-Impregnated	Gingi-Plain #2 Max Soft Twist	none
Retrax 07-920 Size 10/2 Non-Impregnated	Gingi-Plain #3 Max Soft Twist	none